

A NEW GENERATION OF POWER WHEELCHAIRS & SCOOTERS

WHY CHOOSE ELOFLEX?

With an Eloflex, you can regain control of your life and rediscover the freedom of enjoying walks, trips and excursions in your own car. At last, there's a practical and useful power wheelchair, designed in Sweden by wheelchair users, to make your life easier and more comfortable. Our compact electric wheelchairs and hybrid electric scooters can easily be stored in a closet or wardrobe. They let you to use regular taxi services, book flights and trains and continue to use your own car or a rental car while on holiday. Eloflex is exactly what all power wheelchairs should be. Experience the difference with Eloflex and reclaim your freedom!

FOLDABLE

Eloflex is a completely new generation of foldable power wheelchairs with a host of unique features. The modern design lets you fold it, one-handed, into a compact size for handy storage at home or in a car. It takes up little space, and the exceptionally low weight makes it easy to take with you. With Eloflex, you reclaim your freedom – a feeling you can take with you anywhere.



EXCEPTIONAL LOW WEIGHT

Sleek and compact, it's also easy to stow in your car for excursions. At last, you don't have to depend on Paratransit and Disability Transportation to arrive. Eloflex gives you your freedom back.



Eloflex is designed, engineered and manufactured using state-of-theart technology. It's highly efficient Lithium-Ion batteries are far superior to old lead batteries, both in terms of performance and weight. The powerful brushless and gearless direct drive motors are mounted directly on each wheel, making them both quieter and more efficient than the bulky motors used in traditional electric wheelchairs.



ON THE ROAD

Eloflex foldable electric wheelchairs are perfect for trips, excursion by car and taxi rides to your holiday destination. Now you can visit the grandchildren or family and friends on your own terms.

A NEW GENERATION OF POWER WHEELCHAIRS & SCOOTERS

Eloflex is power wheelchair available in models for all types of people. Regardless of your age or physical condition, Eloflex provides added freedom at home and on the go, eliminating the need for external assistance. Its lightweight design makes it easy to pack up and take in your car. At last, you will never have to depend on Paratransit or Disability Transit.

INDEPENDENT FREEDOM

It's easy to understand how the Eloflex solution positively impacts the daily lives of so many people. To start with, it offers a foldable lightweight electric solution that's perfect to take in the car and on trips. More importantly, it offers a new level of mobility that previously did not exist due to the non-foldable nature of traditional power wheelchairs, which often weigh between 80 and 150 kg.

With its compact design and high driving comfort, you have virtually no limits on getting around. You can cruise around the park, or head for town to navigate busy shops, cafes and restaurants.

Modern flight-approved batteries mean that you can check in your Eloflex on planes as well.



COMPACT & VERSATILE

Traditional power wheelchairs are typically large, bulky and weigh more than 100kg. This makes them impactable to fit, let alone lift, into a car for excursions. At home, they take up excessive space and are quite clumsy to drive. But it doesn't have to be that way. With their Eloflex represents a completely modern approach to electric wheelchairs, constructed from aluminum or carbon fiber. With a simple hand grip, it can be effortlessly folded, minimizing the space it occupies and leveraging its exceptionally low weight for easy transport in your car.

Due to its compact dimensions, rear-wheel drive and having the same seat design as a manual wheelchair, Eloflex gives you comparable performance. The tight turning radius and high precision let you to effortlessly maneuver into tight spaces and through narrow doors. Eloflex represents a new generation of power wheelchairs, boasting a lightweight frame and a slim width.







MODEL F

Eloflex F is the market's most sought-after power wheelchair. With unique features like its foldability, compact design, and a lightweight aluminum frame, it has revolutionized the assistive technology sector. The Eloflex Model F has restored freedom to thousands of users who have difficulties walking.

This versatile model is suitable for individuals of all ages, needs and purposes. Its modern technology, featuring powerful brushless motors and dual lithium batteries as standard, provides exceptional driving performance and an impressive range exceeding 30 km.



 WEIGHT IN KG
 RANGE IN KM
 SEAT HEIGHT IN CM
 MAX LOAD IN KG
 EFFECT IN WATT

 26
 30
 52
 120
 2X 250

MODEL L

Eloflex L is a compact, no-nonsense power wheelchair for "no fuss" people on the move. It conveniently folds in seconds, making it easy to take with you when traveling by car or by plane.

The Eloflex Model L is an excellent option for those who want ease and peace of mind. With solid wheels, there's no need to inflate them and you can say goodbye to concerns about flat tires. Its sleek and nimble format is particularly advantageous in indoor settings, allowing you to effortlessly navigate through narrow doorways and confined spaces.



WEIGHT IN KG RANGE IN KM SEAT HEIGHT IN CM MAX LOAD IN KG EFFECT IN WATT

25,5 30 52 120 2X 200

MODEL P

Eloflex P has a built-in electrical tilt function in the seat, designed to facilitate entry and exit from the chair. This model is intended for people who face difficulties in standing and getting out of a standard chair. With the electric lift, the transition is both easier and safer. Not only does the electric seat lift makes the transition easier but also safer. And it simplifies reaching higher places, such as cabinets and shelves, both at home and in stores.

Model P features a higher seat height compared to our other models, further easing the exit process. As a result, Eloflex P is an excellent choice for slightly taller users, offering a higher seat position and the flexibility to raise the backrest seamlessly.



 WEIGHT IN KG
 RANGE IN KM
 SEAT HEIGHT IN CM
 MAX LOAD IN KG
 EFFECT IN WATT

 28
 30
 54
 160
 2X 250

MODEL D2

Eloflex D2 is an innovative foldable power wheelchair, sporting the same sleek design as our other power wheelchairs but with a smart twist – it can be disassembled into two parts. These two parts provide a loading weight of 14 and 16 kg, making it easier to transport and load the chair into your car trunk.

The stable D2 model accommodates users weighing up to 160 kg and offers comfort and ease of maneuverability. Eloflex D2 empowers you to live an active life, whether it's going for walks, traveling or taking everyday excursions.



160

2X 250

WEIGHT IN KG RANGE IN KM SEAT HEIGHT IN CM MAX LOAD IN KG EFFECT IN WATT

14 + 16 30 50

MODEL H

With Eloflex Model H, we have done the impossible – it's a powerful electric wheelchair specially designed for our heavier and taller users, capable of accommodating users weighing up to 200 kg. Remarkably, it can also be folded and transported in a regular car.

Model H is a spacious chair with great driving characteristics and plenty of room for fully grown individuals. This model offers flexible use thanks to its generous seat width, comfortable seat height, and a design that provides a comfortable driving position, even for taller individuals.



 WEIGHT IN KG
 RANGE IN KM
 SEAT HEIGHT IN CM
 MAX LOAD IN KG
 EFFECT IN WATT

 43
 40
 54
 200
 2X 400

MODEL K

Eloflex K is a compact, lightweight, and foldable power wheelchair designed for children and teenagers aged 7 to 15. Model K empowers children and young adults to actively engage in their surroundings, making it easier for them to participate in everyday activities such as school trips, shopping, dining, excursions and visiting friends and family – all without the need for specialized transportation.

The K model boasts a compact and lightweight design, featuring the same elegant aesthetics found in other Eloflex models. The chair is thoughtfully crafted to adapt as its users grow, offering simple height and width adjustments to accommodate a child's evolving requirements. Model K is straightforward to drive, making it possible for children as young as 7 to navigate independently.



WEIGHT IN KG RANGE IN KM SEAT HEIGHT IN CM MAX LOAD IN KG EFFECT IN WATT

26 30 47-52 90 2X 250

MODEL Z

Eloflex Z stands as the most exclusive foldable power wheelchair available, both in our own range and on the market in general. This model is tailored for those who want the "best of the best" and hold their assistive device to the highest standards. Crafted from carbon fiber, the frame ensures an exceptionally low weight of only 23kg. Its format is compact, flexible, and foldable.

The meticulously design features and ingenious details offer the most advanced adjustment options on the market, allowing for optimal individualization. Unrivaled driving performance, an extraordinary range of up to 45 km, and exceptional reliability are just a few examples of the model's unique qualities. Eloflex Z is the electric wheelchair of choice for those who insist on the best and refuse to compromise on anything less than excellence.



WEIGHT IN KG	RANGE IN KM	SEAT HEIGHT IN CM	MAX LOAD IN KG	EFFECT IN WATT
			• • • • • • • • • • • • • • • • • • • •	
20	30	52	120	2X 250

MODEL C3

The Eloflex C3 combines the most cherished features of an Eloflex with the cutting-edge advantage of carbon fiber construction, making it incredibly lightweight and your perfect companion for daily use.

Weighing just 20 kg, the Eloflex C3 boasts a sleek design complemented by sturdy wheels and powerful Eloflex motors, ensuring a superior and smooth riding experience. With a maximum user weight capacity of 160 kg and additional space for extra batteries, the Eloflex C3 is engineered to meet and surpass all your expectations, offering both reliability and enhanced mobility.



WEIGHT IN KG RANGE IN KM SEAT HEIGHT IN CM MAX LOAD IN KG EFFECT IN WATT

20 **30** 55 160 2X 250

MODEL X

Eloflex X is a unique design tailored specially for our elderly users. It stands as the most popular senior power wheelchair on the market, built on the well-proven and widely appreciated crossbrace frame technology. Folding the chair is a breeze – simply lift the seat.

The Eloflex X is the result of years of invaluable feedback, comments and requests received from our loyal elderly users. It represents a "pure senior model" with all the functions and solutions one could desire.



 WEIGHT IN KG
 RANGE IN KM
 SEAT HEIGHT IN CM
 MAX LOAD IN KG
 EFFECT IN WATT

 25
 25
 51
 120
 2X 250

MODEL S1

With Eloflex's groundbreaking and oneof-a-kind hybrid scooter, you can move seamlessly both indoors and outdoors, all while easily transporting your scooter in your car or taking it on a flight. So there's nothing stopping you to take long walks, shop, meet friends and family, and travel.

A lightweight scooter that's easy to bring along, yet robust enough for outdoor adventures in challenging environments. The Eloflex S1 is built to withstand all seasons, including the harsh Nordic winters, making it a versatile and reliable choice for year-round mobility.



WEIGHT IN KG RANGE IN KM SEAT HEIGHT IN CM MAX LOAD IN KG EFFECT IN WATT

37,5 30 55-60 120 2X 250

ACCESSORIES

Make your Eloflex even more brilliant by customizing it to your individual needs with our popular and smart accessories.



ELOFLEX PREMIUM CUSHIONS

Ease by Eloflex offers a range of Swedish-made ergonomic cushions designed to optimize your Eloflex user experience.



SIDE BAG

With the Eloflex Side Bag, you can conveniently store items like your wallet, phone, and other valuables for easy access.



PROTECTIVE TRAVEL COVER

Our Protective Travel Cover is the perfect accessory for Eloflex power wheelchairs enthusiasts who are fussy about their chairs and want to take extra precautions when transporting them.



ADJUSTABLE HEADREST

The Eloflex Headrest is a highly versatile option for extra head and neck support, designed to increase comfort and usability for users with different needs.



ATTENDANT CONTROL SYSTEM

Thanks to our Attendant Control System, individuals walking behind or beside the wheelchair can drive the electric wheelchair with the help of the joystick, allowing the person seated to effortlessly ride along.



JOYSTICK EXTENSIONS

Eloflex offers two variants of joystick extensions to suit a wide range of needs.



CHARGING ADAPTER

The Eloflex Charging Adapter allows you to directly connect the charger to the battery.



ADJUSTABLE FOOTREST

Our Adjustable Footrest allows you to fine-tune the angle of the foot plate anywhere from 0 to 90 degrees.



ELOFLEX RAMPS

Easily load your Eloflex into the car with the help of load ramps that are collapsible and easy to lift.

COMFORT FOR EVERY OCCASION WITH EASE BY ELOFLEX

Prolonged periods of sitting can sometimes lead to discomfort in the back and buttocks, skin irritation, cramps or, in the worst case, pressure sores and ulcers. The ergonomic products in our Ease range help prevent unnecessary discomfort and avoid pressure ulcers. They also provide relief from pain

and alleviate existing issues. With a variety of products to choose from, you can easily get the best support for your arms, neck and back whether at the kitchen table, in front of the TV, in the car, or in a manual wheelchair.



PREMIUM CUSHION - EASE ONE

Ease One is a soft and ultra-comfortable premium cushion designed for people who spend extended periods sitting. The cushion features a bed of coil springs strategically placed under the buttock bones and the tailbone. The contours are also shaped to enhance support and stability for the thighs. Ease One is particularly well-suited for use in electric wheelchairs as its advanced suspension protects the back and tailbone on uneven surfaces.

PREMIUM CUSHION - EASE WEDGE

Ease Wedge is a wedge-shaped seat cushion that is raised at the front. It is specifically designed to enhance the comfort of individuals with longer legs who are seated for longer periods of time. The ergonomic design distributes pressure and offers a soft feeling that relieves both the femurs, the pelvic "sitting bones," and the tailbone.

EASE ARMREST CUSHION

Ease Armrest Cushion is a plush armrest cover designed for your Eloflex. Featuring 1 cm thick, soft foam padding, it protects you from vibrations while reducing the risk of pressure sores and irritation. It also improves stability, preventing your arm from slipping off the armrest.

EASE BACKREST TUBE PADDING

Ease Backrest Tube Padding is a soft foam covering designed for the backrest of your Eloflex. It offers a softer and more pleasant feeling when you lean back in your chair. You also get extra back support and comfort when powering along on uneven surfaces.

EASE LUMBAR SUPPORT

Ease Lumbar Support is designed in firm and durable viscoelastic "memory" foam. Unlike regular foam, memory foam reacts to body heat and pressure, becoming softer and contouring to the lumbar "lower" spine to reduce back fatigue.

EASE INCONTINENCE COVER

Dealing with incontinence while seated in a wheelchair shouldn't be stressful or complicated. Our waterproof incontinence cover simplifies and safeguards your everyday life, providing excellent protection against rain and moisture when you're outdoors.

MARIE FOUND TRUE FREEDOM WITH ELOFLEX



Hopptimisten, or Marie Hav Lunkvist as she is actually called, has a background as an active businesswoman who is currently in a wheelchair due to illness. She sees life as a gift and is filled with both hope and optimism!

I got my Eloflex in February 2018, I had then had other wheelchairs for about six years. I am an active person who travels a lot with lectures and I'm involved in many different community assignments both privately and with the family. I travel to both the mountains and go on archipelago trips, and it became absolutely necessary for me to have a flexible and lightweight electric wheelchair on these occasions.

It is not only when I travel that the Eloflex is with me, but everything from not getting into crowded shops or such an everyday thing as the neighbour's door being too narrow for my regular electric wheelchair, is a memory now. I no longer have to worry about not being able to do what I want, when I want to, and on my own terms. That's freedom for real!



TECHNICAL SPECIFICATION

MODEL	F	L	P	D2	н	K	С	Z	C3	X	S 1
NET WEIGHT	27 kg	25,5 kg	30 kg	14+16 kg	42,5 kg	28,5 kg	16 kg	23 kg	20 kg	25 kg	37,5 kg
RANGE	30 km	30 km	30 km	30 km	40 km	30 km	15 km	30 km	30 km	25 km	30 km
MAX LOAD	120 kg	120 kg	160 kg	160 kg	200 kg	90 kg	100 kg	120 kg	160 kg	120 kg	120 kg
SPEED	6,5 km/h	6,5 km/h	6,5 km/h	6,5 km/h	8 km/h	6,5 km/h	6 km/h	6,5 km/h	6,5 km/h	6,5 km/h	9 km/h
BATTERY	2x Litium ION 10AH/24V	2x Litium ION 10AH/24V	2x Litium ION 10AH/24V	2x Litium ION 10AH/24V	2x Litium ION 20AH/24V	2x Litium ION 10AH/24V	2x Litium ION 12AH/24V	2x Litium ION 10AH/24V	2x Litium ION 10AH/24V	2x Litium ION 10AH/24V	2x Litium ION 10AH/24V
TURNING RADIATION	78 cm	76 cm	85 cm	89 cm	99 cm	78 cm	60 cm	78 cm	87 cm	88 cm	100 cm
WHEELS FRONT	8" solid	8" solid	8" solid	8" solid	9" air	8" comfort	7" solid	8" comfort	10" air	8" solid	12" air
WHEELS REAR	12" air	10" solid	12″ air	12" air	13″ air	12" air	8,5" solid	12" air	12" air	12" air	12″ air
SEAT HEIGHT	52 cm	52 cm	54 cm	50 cm	54 cm	47-52 cm	47,5 cm	52 cm	55 cm	51 cm	55-60 cm
SEAT WIDTH	45/50/55 cm	45/50/55 cm	45/50/55 cm	45/50/55 cm	51/56/61 cm	36-44 cm	44 cm	49 cm	43 cm	45/50/55 cm	45/50/55 cm
SEAT DEPTH	43-48 cm	43-48 cm	43-48 cm	44-47 cm	51 cm	36 cm	42 cm	45 cm	44 cm	43-48 cm	40 cm
SIZE (HXWXD)	95x63x100 cm	95x59x99 cm	90x62x107 cm	96x62x107 cm	103x66x 115 cm	89x63x 100 cm	93x57x 93 cm	100x62x 110 cm	95x63x100 cm	100x62x 102 cm	91x63x 105 cm
SIZE FOLDED (HXWXD)	75x63x40 cm	70x59x35 cm	74x62x40 cm	67/52x 63/63x 50/35 cm	83x66x45 cm	75x55x36 cm	79x57,5x 32,5 cm	75x64x37 cm	75x63x40 cm	76x32x78 cm	55x63x105 cm

Net weight don't include batteries. Recommended air pressure for air-filled wheels = 2.5 kg. Eloflex is classified as a medical device with the Swedish Medicines Agency and meets the applicable requirements in the EU regulation (2017/745), Medical Device Regulations (MDR), which ensures the safety and performance of medical devices. It is CE approved, tested and approved according to EN12184. ISO-7176. IEC-60601 and the batteries are flight approved according to UN38.3.

Visit www.eloflex.eu for more information, videos, pictures and your nearest reseller.



Eloflex HQ

Hammarbacken 4B 191 49 Sollentuna SWEDEN